

FEMORAL

Indications for Use: The Medcomp® Femoral Catheter is designed for acute hemodialysis and apheresis. It may be inserted percutaneously and the preferred site is the femoral vein. This catheter is indicated for a duration less than (30) days.

Contraindications: This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane

Material

Tapered Tip

Tip Design

7F, 8.5F, 9F

French Sizes

**Single Lumen,
Straight**

Configurations



MULTIPLE CONFIGURATIONS
Available in double or single lumen designs.

**USEFUL IN OVER-THE-WIRE AND
OVER-THE-NEEDLE INSERTIONS**

CATHETER ONLY		
MCF64	7F X 13.5CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCF55	8.5F X 14CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCF10	8.5F X 25CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCF6	8.5F X 15CM FEMORAL CATHETER, STRAIGHT	10/BOX

BASIC SET		
MCFK64	7F X 13.5CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCFK100	8.5F X 25CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCFK550	9F X 14CM FEMORAL CATHETER, STRAIGHT	10/BOX
MCFK60	8.5F X 15CM FEMORAL CATHETER, STRAIGHT	10/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) 2 End Cap(s)
- (1) "Y" Hub

*MCF6 Only

BASIC SET CONTAINS:

- (1) Catheter
- (1) 2 End Cap(s)
- (1) "Y" Hub
- (1) Introducer Needle
- (1) Guidewire

*MCFK6 Only

[Refer to the Table of Contents for Important Risk Information regarding this device.](#)