



TRIO-CT®

Indications for Use: The Trio-CT® Triple Lumen Catheter is indicated for use in attaining short-term (less than 30 days) vascular access for hemodialysis and apheresis. The third internal lumen is intended for infusion, power injection of contrast media and central venous pressure monitoring. The catheter is intended to be inserted in the jugular, femoral or subclavian vein as required. The maximum recommended infusion rate is 5ml/sec for power injection of contrast media.

Contraindications: This catheter is intended for short-term (less than 30 days) vascular access only and should not be used for any purpose other than indicated in these instructions. This device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient to accommodate the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors may prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane	Tapered Tip	13.5F	Straight
Material	Tip Design	French Size	Configuration



TRIPLE INTERNAL
LUMEN DESIGN



CATHETER ONLY		
AS01135120E.	13.5F X 12CM TRIO-CT® CATHETER ONLY	10/BOX
AS01135150E.	13.5F X 15CM TRIO-CT® CATHETER ONLY	10/BOX
AS01135200E.	13.5F X 20CM TRIO-CT® CATHETER ONLY	10/BOX
AS01135240E.	13.5F X 24CM TRIO-CT® CATHETER ONLY	10/BOX
AS01135300E.	13.5F X 30CM TRIO-CT® CATHETER ONLY	10/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Dilator
- (1) Attachable Suture Wing
- (3) End Caps
- (1) Patient ID Card
- (1) Patient Chart Sticker

BASIC SET		
AS01135121E.	13.5F X 12CM TRIO-CT® BASIC SET	5/BOX
AS01135151E.	13.5F X 15CM TRIO-CT® BASIC SET	5/BOX
AS01135201E.	13.5F X 20CM TRIO-CT® BASIC SET	5/BOX
AS01135241E.	13.5F X 24CM TRIO-CT® BASIC SET	5/BOX
AS01135301E.	13.5F X 30CM TRIO-CT® BASIC SET	5/BOX

BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Attachable Suture Wing
- (3) End Caps
- (1) Patient ID Card
- (1) Patient Chart Sticker

BASIC SET ME		
AS01135121ME	13.5F X 12CM TRIO-CT® BASIC SET	5/BOX
AS01135151ME	13.5F X 15CM TRIO-CT® BASIC SET	5/BOX
AS01135201ME	13.5F X 20CM TRIO-CT® BASIC SET	5/BOX
AS01135241ME	13.5F X 24CM TRIO-CT® BASIC SET	5/BOX
AS01135301ME	13.5F X 30CM TRIO-CT® BASIC SET	5/BOX

BASIC SET ME CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Scalpel
- (1) Adhesive Wound Dressing
- (1) Attachable Suture Wing
- (3) End Caps
- (1) Patient ID Card
- (1) Patient Chart Sticker

Refer to the Table of Contents for Important Risk Information regarding this device.