



TRIO-CT®

Indications for Use: The Trio-CT® Triple Lumen Catheter is indicated for use in attaining short-term (less than 30 days) vascular access for hemodialysis and apheresis. The third internal lumen is intended for infusion, power injection of contrast media and central venous pressure monitoring. The catheter is intended to be inserted in the jugular, femoral or subclavian vein as required. The maximum recommended infusion rate is 5ml/sec for power injection of contrast media.

Contraindications: This catheter is intended for short-term (less than 30 days) vascular access only and should not be used for any purpose other than indicated in these instructions. This device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient to accommodate the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors may prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane	Tapered Tip	13.5F	Straight
Material	Tip Design	French Sizes	Configuration



TRIPLE ROUND INTERNAL LUMEN DESIGN



CATHETER ONLY		
AS01135120	13.5F X 12CM TRIO-CT®	10/BOX
AS01135150	13.5F X 15CM TRIO-CT®	10/BOX
AS01135200	13.5F X 20CM TRIO-CT®	10/BOX
AS01135240	13.5F X 24CM TRIO-CT®	10/BOX
AS01135300	13.5F X 30CM TRIO-CT®	10/BOX

CATHETER ONLY CONTENTS:

- (1) Catheter
- (2) Dilators
- (1) Attachable Suture Wing
- (3) End Caps
- (1) Patient ID Card
- (1) Chart Sticker

BASIC SET		
AS01135121	13.5F X 12CM TRIO-CT®	5/BOX
AS01135151	13.5F X 15CM TRIO-CT®	5/BOX
AS01135201	13.5F X 20CM TRIO-CT®	5/BOX
AS01135241	13.5F X 24CM TRIO-CT®	5/BOX
AS01135301	13.5F X 30CM TRIO-CT®	5/BOX

BASIC SET CONTENTS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Adhesive Wound Dressing
- (1) Attachable Suture Wing
- (3) End Caps
- (1) Patient ID Card
- (1) Chart Sticker

MAX BARRIER TRAY		
AS01135123	13.5F X 12CM TRIO-CT®	3/BOX
AS01135153	13.5F X 15CM TRIO-CT®	3/BOX
AS01135203	13.5F X 20CM TRIO-CT®	3/BOX
AS01135243	13.5F X 24CM TRIO-CT®	3/BOX
AS01135303	13.5F X 30CM TRIO-CT®	3/BOX

MAX BARRIER TRAY CONTENTS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Adhesive Wound Dressing
- (1) Attachable Suture Wing
- (3) End Caps
- (3) 4x4 Gauze
- (2) Chloraprep
- (3) 10cc Syringe
- (1) Filter Straw
- (1) Lidocaine/Saline Sticker Sheet
- (1) 3-0 Silk Curved Suture
- (2) Safety Hypodermic Needles
- (1) Hemostat
- (1) Full Body Max Barrier Drape
- (1) XL Gown
- (1) Bouffant Cap
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Steri-Strips
- (1) Face Mask
- (1) Surgical Face Mask
- (1) Patient ID Card
- (1) Chart Sticker

Refer to the Table of Contents for Important Risk Information regarding this device.